

BUNSEN BURNER STRIP STORY

The sentences below are muddled up. Put them in the correct order to find the instructions for using a Bunsen burner safely.

- 1.
- 2.
- 3.
- 4.
- 5.
- 6.
- 7.
- 8.
- 9.
- 10.

Have a lighted splint ready.

Check that the hose fits tightly or it could come off.

The Bunsen burner then goes on top of the board.

USING A BUNSEN BURNER SAFELY

Turn the gas tap on full.

Attach the rubber hose of the Bunsen to the gas tap.

Leave the collar closed until you use the Bunsen burner.

Close the air hole of the Bunsen before lighting.

First you must put a fireproof board down on the bench.

Light the Bunsen about five centimetres above the barrel.