

TAJUK : ISIPADU

NAMA :

TONTON VIDEO DI BAWAH :



B. Tolak ukuran isi padu dengan mengumpul semula



Cari beza isi padu susu kotak dalam unit ml.

$$1.1 \ell - 240 \text{ ml} = \text{ } \text{ml}$$

$$\begin{array}{r}
 \begin{array}{l} 1.1 \ell \\ = 1 \ell 100 \text{ ml} \end{array} \rightarrow \begin{array}{r} \ell \quad \text{ml} \\ 1 \quad 00 \\ - \quad 240 \\ \hline \end{array} \leftarrow \begin{array}{r} 10 \quad 10 \\ + \quad + \quad 0 \quad 0 \text{ ml} \\ - \quad 240 \text{ ml} \\ \hline \end{array}
 \end{array}$$

$$1.1 \ell - 240 \text{ ml} = \text{ } \text{ml}$$

Isi padu kolam
= 101 ℓ



Isi padu air yang ada
= 89 ℓ 800 ml

Kira isi padu air yang diperlukan untuk memenuhi kolam.

$$101 \ell - 89 \ell 800 \text{ ml} = \text{ } \ell \text{ } \text{ ml}$$

$$\begin{array}{r} \ell \quad \text{ml} \\ 101 \quad 000 \\ - 89 \quad 800 \\ \hline \end{array}$$

$$101 \ell - 89 \ell 800 \text{ ml} = \text{ } \text{ ml}$$

Rajah menunjukkan jumlah isi padu air tebu di dalam 2 buah jag.
1 ℓ 890 ml air tebu telah dituang ke dalam kole.
Hitung isi padu air tebu yang tinggal dalam unit ml.



2 ℓ 830 ml

$$2 \ell 830 \text{ ml} - 1 \ell 890 \text{ ml} = \text{ } \text{ ml}$$

$$\begin{array}{r} \ell \quad \text{ml} \\ 2 \quad 830 \\ - 1 \quad 890 \\ \hline \end{array}$$

$$\begin{array}{r} 17 \quad 13 \\ + 8 \quad 30 \text{ ml} \\ - 8 \quad 90 \text{ ml} \\ \hline \text{ } \text{ ml} \end{array}$$

$$2 \ell 830 \text{ ml} - 1 \ell 890 \text{ ml} = \text{ } \text{ ml}$$

Jadual menunjukkan isi padu susu di dalam sebuah bekas sebelum dan selepas sebahagiannya digunakan untuk membuat minuman berkhasiat.

	Sebelum	Selepas
Isi padu susu	2 ℓ 40 ml	750 ml

Berapakah isi padu susu yang digunakan dalam unit ℓ dan ml?

$$2 \text{ ℓ } 40 \text{ ml} - 750 \text{ ml} = \text{ } \text{ ℓ } \text{ } \text{ ml}$$

$$\begin{array}{r} \text{ℓ} \quad \text{ml} \\ 2 \quad 40 \\ - \quad 750 \\ \hline \end{array}$$

$$\begin{array}{r} 2 \text{ ℓ } 40 \text{ ml} \\ - 750 \text{ ml} \\ \hline \end{array}$$



$$2 \text{ ℓ } 40 \text{ ml} - 750 \text{ ml} = \text{ } \text{ ml}$$