

Boletín.- Cambio de Unidades

1) Completa

a. $600\text{dm} = \underline{\hspace{2cm}}\text{hm}$

d. $400\text{cm} = \underline{\hspace{2cm}}\text{dam}$

b. $410\text{dam} = \underline{\hspace{2cm}}\text{km}$

e. $0,009\text{km} = \underline{\hspace{2cm}}\text{m}$

c. $0,34\text{dm} = \underline{\hspace{2cm}}\text{mm}$

f. $200\text{hm} = \underline{\hspace{2cm}}\text{dam}$

2) Completa

a. $0,05\text{kl} = \underline{\hspace{2cm}}\text{dl}$

d. $1000\text{dal} = \underline{\hspace{2cm}}\text{hl}$

b. $340\text{ml} = \underline{\hspace{2cm}}\text{cl}$

e. $75\text{cl} = \underline{\hspace{2cm}}\text{l}$

c. $0,9\text{l} = \underline{\hspace{2cm}}\text{dl}$

f. $300\text{ml} = \underline{\hspace{2cm}}\text{cl}$

3) Completa

a. $80\text{dg} = \underline{\hspace{2cm}}\text{g}$

d. $700\text{dg} = \underline{\hspace{2cm}}\text{hg}$

b. $150\text{mg} = \underline{\hspace{2cm}}\text{dg}$

e. $4000\text{mg} = \underline{\hspace{2cm}}\text{kg}$

c. $0,03\text{dg} = \underline{\hspace{2cm}}\text{dag}$

f. $0,9\text{dg} = \underline{\hspace{2cm}}\text{cg}$

4) Completa

a. $2.422\text{m} = \underline{\hspace{2cm}}\text{km}$

d. $400\text{mm} = \underline{\hspace{2cm}}\text{dm}$

b. $4,2\text{dam} = \underline{\hspace{2cm}}\text{hm}$

e. $1,6\text{km} = \underline{\hspace{2cm}}\text{dam}$

c. $256\text{cm} = \underline{\hspace{2cm}}\text{m}$

f. $0,52\text{km} = \underline{\hspace{2cm}}\text{m}$

5) Completa

a. $1,56\text{kl} = \underline{\hspace{2cm}}\text{dal}$

d. $324\text{l} = \underline{\hspace{2cm}}\text{hl}$

b. $5,4\text{hl} = \underline{\hspace{2cm}}\text{dl}$

e. $75\text{cl} = \underline{\hspace{2cm}}\text{l}$

c. $0,58\text{l} = \underline{\hspace{2cm}}\text{dl}$

f. $234\text{ml} = \underline{\hspace{2cm}}\text{dl}$

6) Completa

a. $3\text{kg} = \underline{\hspace{2cm}}\text{g}$

d. $28,7\text{dg} = \underline{\hspace{2cm}}\text{g}$

b. $1,4\text{hg} = \underline{\hspace{2cm}}\text{dag}$

e. $470\text{mg} = \underline{\hspace{2cm}}\text{cg}$

c. $420\text{g} = \underline{\hspace{2cm}}\text{kg}$

f. $39\text{dg} = \underline{\hspace{2cm}}\text{mg}$